
50 50 T

50 50 T

Key points. *Hauptmerkmale.* Puntos clave. *Points clé.* Punti chiavi. Pontos chave.

1 Connector: MP-5.
Anschluss: MP-5.
Conector: MP-5.
Connecteur: MP-5.
Connettore: MP-5.
Conector: MP - 5.

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pH electrode 50 50 T

Introduction

The **50 50 T** is a universal pH electrode, with built-in A.T.C. Pt 1000. These are low maintenance electrodes, as they do not require monitoring of electrolyte levels.

SAFETY NOTE

Before use with hazardous substances, find out about the relevant hazards, protective measures, and what to do in an emergency.

Electrode preparation

Remove, by unscrewing it, the storage protector containing electrolyte ④.

 *MP-5 connector (IP 67). Connect and disconnect by pushing/pulling. Do not unscrew it.*

Calibration

 *See pH-meter instruction manual.*

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases. In field measurements, these types of electrodes allow stirring the sample directly holding the electrode by the handle.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

 *Do not leave the electrode immersed in distilled water.*

 *Do not bend the cable on the handle level.*

 *Never unscrew the MP-5 connector.*

Maintenance

- Maintain the membrane's surface ⑨ clean.
- When an electrode gets contaminated, the simply rinsing with distilled water is not sufficient. The most suitable cleaning solution is which acts in a most selective way over the contamination (see regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrodes regenerated solutions.

Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode, cable or connector.

Substitute the electrode.

... fortuitous reading.

Broken cable.

Substitute the electrode.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm and/or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or replace the electrode.

Storage

Place the electrode in the storage protector ④ with appropriate electrolyte KCl 3M.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

 *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

pH measuring range	0 to 14
Operating temp.	0 to 80 °C
Reference element	Ag wire coated with AgCl
Diaphragms	ceramic
Electrolyte	gel
Body material	polycarbonate
Min. immersion depth	18 mm 
Temperature sensor	Pt 1000

Applications

Aqueous samples in general.

Limitations

Low-conductivity, dirty or viscous solutions.

Products with colloids or solids in suspension.

Samples containing sulphides, reducing sugars or other substances which react with the silver ions.

Temperatures > 80 °C.

QUALITY CERTIFICATE

The attached electrode 50 50 T has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4 to 7 (at 25 °C) > 98%.
- Response time, pH 4 to 7 < 20 s.
- Temperature error ≤ 0.5 °C.

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